MICROSCOPIOS DE OPERACIONES OPERATION MICROSCOPES

Serie / Series 191

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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WARRANTY

AUXILAB, S.L. warrants this equipment for a period of twelve months from the date of purchasing to be free from any defects of manufacturing that affects its functioning, we will replace these pieces which break due to a manufacturing defect. The change of these pieces will be free during the warranty period and will not extend it.

This warranty does not extend to any products which have been subject to misuse, have been manipulated by somebody else other than our Technical Service, or have not been installed properly. This also affects to equipments that have not used the original accessories.

This limited warranty does not extend to any equipments which have been damaged as a result of catastrophic causes (fire, flood,...), atmospheric, blows or falls.

The costs caused by this repair, the transport of the equipments to our Technical Services and its return will be paid by the user.

It is necessary the presentation of the Warranty Certificate properly filled (without strokes or crossing outs) in order to achieve any free of charge operation of our Technical Service during the warranty period. This warranty does not extend to the maintenance of the equipment or its pieces, or the damages produced by the everyday use of the apparatus.

To provide a better service, when you send an equipment to our Technical Service we would be grateful if you could provide us with the following information.

Supplier:			
Equipment:		Model:	
Serial number:			
Date of purchase:			
Brief description of the damage:			

AUXILAB, S.L.

Polígono Morea Norte, 8 - 31191 BERRAIN (Navarra)

Tel(948) 310 513 - Fax(948) 310 500 - 312 071

e-mail: asistencia@auxilab.com

<http://www.auxilab.es>

8. TROUBLESHOOTING

Check the following points before sending the microscope to our Technical Assistance Department:

If the lamp does not switch on:

- ◆ Check it is properly connected to the right current intake.
- ◆ Check both fuse and lamp are in good conditions. If not you should replace the fuse.

If there is dust or dirt in the field of vision:

Check there is dust on the objectives or upper lens of the eyepiece. Once you find out where the dirt is you should clean it as explained before.

9. ACCESSORIES

There are two accessories available for Zuzi operation microscopes:

Wall support for microscope, code 90190901

- ◆ With this wall support great manoeuvres can be performed as long as the microscope is used at the same place and it is not being moved.

Desktop support for microscope, code 90190902


- ◆ It enables the microscope use in places where it has not to be moved.

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6. USEFUL ADVICE

You must clean lenses after using the microscope. **Never use a chamois duster!** Cleaning has to be done with a special paper for cleaning glasses or a linen cloth. You should never immerse objectives in any liquid, as it may provoke lens detachment.

Possible causes of a deficient image

In case you see a deficient image you should check the following:

- ◆ Illumination is well done and luminous intensity is neither excessive nor too weak.
- ◆ Cleanliness of all optical system, proceeding as follows: rotate the eyepieces and in case dust specks move, clean them.

7. MAINTENANCE

- ◆ All the processes of use mentioned below will not have any value unless you keep a continued and careful maintenance.
- ◆ The essential maintenance of these microscopes is the same for any optical equipment. They should always be kept clean and covered to avoid dust.
- ◆ Lenses must not be disassembled by the user. Were there any dust or dirt to be cleaned, you should clean it with a smooth piece of cloth, fluff-free.
- ◆ To remove dust from lenses you should blow them with a plastic bulb or clean them with a natural horsehair brush.
- ◆ When cleaning mechanical or plastic parts you should not use abrasive solvents that could damage either paint or finish. You should always use neutral detergent.
- ◆ In case the lamp blows you should replace it for another one, making sure it is a 12V 100W Zuzi original one. Please take care and do not touch the lamp with bare hands.

IMPORTANT: Before replacing the lamp make sure you have disconnected the microscope from the electric intake. Check you are using 12V 100W Zuzi original halogen lamps, as other types can provoke malfunction.

- ◆ To change the lamp, press the two black buttons next to the lighting wire input (1.17), then pull rotating the bulb a little. To replace it you should have to wait until it cools down, move the spring that tightens the bulb and take it out. When pulling the lamp take out the lamp's arms from the lamp holder and replace it.
- ◆ To change the fuse unscrew the *Fuse* button (1.16), and then replace it.

4. SETTING UP

- ◆ Take out of the package the stand and all the rest of pieces of the microscope and check it does not present any damage because of the shipment.
- ◆ Put the microscope stand on top of a flat, horizontal surface and keep it away from dust, humidity, heat and vibrations.
- ◆ Assemble the different pieces of the microscope according to the picture. Fix the base with the brakes on the wheels (1.10), place the column with the lower screw, then assemble the piece with the electric components and fix it with the screws (1.13) and (1.18). Now assemble the extending arm and the head on it, and put the video adapter if needed. There is a photographic attachment as an accessory.

5. INSTRUCTIONS OF USE

- ◆ Plug the microscope to the net -mind that the switch located on the lower part (1.12) is on the proper position-. Then, switch on the microscope by pressing the button (1.15) located on the upper part of (1.12).
- ◆ Regulate height and position of the extending arm according to the height of the sample to be observed.
- ◆ Adjust head position by using its handlebars (1.5). You can move the head to the desired position thanks to the movable arm.
- ◆ Personalize the interpupillary distance by adjusting the eyepieces to the desired distance. You will know it when there is a total fusion of the two images.
- ◆ Notice that during the observation you will need to make little adjustments such as dioptre compensation on eyepiece tubes. By closing alternatively one eye, then the other, you will appreciate a certain focusing difference. In order to correct this fault you should move dioptre compensation rings located on the eyepieces. To do so, you should first focus using only one eyepiece and setting its dioptre adjust to zero, then graduate the other eyepiece using its dioptre compensation ring.

Note: Note: It is advisable to memorize these values once you find them out, most of all if the microscope is shared by more than one user, so as to avoid repeating this process every time the microscope is being used.

- ◆ The focusing is done by using the focusing knobs (1.7). There are four different focusing distances depending on the space needed. Change the objective according to your needs, working distances are marked on the objectives (f').

Note: Fixing screws, focusing knobs and handlebars are sterilizable to avoid contamination.

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Mechanical specifications

STAND:

- ♦ Stand with 5 wheels, very stable. Two of the wheels have a built in brake to provide stability to the microscope.
- ♦ Equipped with an extending and rotary arm.

HEAD:

- ♦ Eyepiece holders are provided with dioptre adjustment and interpupillary distance adjustment.

♦ LIGHTING:

Halogen lamp for cold coaxial lighting 12V/ 100W, built in power source (220-230V 50Hz.).

Note: The microscope has the power point either at 110V or 220V, so please check before connecting the equipment that the switch, located on the lower side of the microscope body, is on the proper position.

Optical specifications

Description	VALUES											
	Objective I			Objective II			Objective III			Objective IV		
	f' = 200 mm			f' = 250 mm			f' = 300 mm			f' = 400 mm		
Magnification	1.6x	1x	0.63x	1.6x	1x	0.63x	1.6x	1x	0.63x	1.6x	1x	0.63x
Total magnif.	12x	7.5x	4.7x	9.6x	6x	3.7x	8x	5x	3x	6x	3.8x	2.3x
Field of view*	15.8	25.3	40.6	19.7	31.6	50.7	23.6	37.9	60.8	31.5	50.5	81
Light* exit diameter	1.04	1.66		1.04	1.66		1.04	1.66		1.04	1.66	
Light exit distance	15.3 mm											
Resolution (LP/mm)	67	44.5	29.7	60	35.4	27	47.2	31.5	23.6	35	27	19.8

f' = Working distance

* Distance in mm


1. USES OF THE INSTRUMENT

Zuzi operation microscopes are frequently used in Microsurgery, especially in fields such as Orthopaedic Surgery, Cardiovascular Surgery, suture of blood vessels, etc.

These microscopes are widely used in veterinary as well as in research and are equipped with a magnifying system which provides three different magnifications. With adjustment control and four interchangeable objectives so as to work with various working distances.

2. DESCRIPTION

- 1.1 Head with dioptre compensation ring and interpupillary distance adjustment
- 1.2 Video and photographic attachment
- 1.3 Photographic attachment (optional as an accessory)
- 1.4 Adjustment screw
- 1.5 Handlebars to move the head up and down
- 1.6 Filters
- 1.7 Focusing knobs
- 1.8 Objectives
- 1.9 Stand with wheels
- 1.10 Non-sliding system
- 1.11 Support column
- 1.12 Electronic component
- 1.13 Fixing screw
- 1.14 Light intensity control
- 1.15 On/off switch
- 1.16 Fuse
- 1.17 Lighting wire input
- 1.18 Fixing screw
- 1.19 Arm tension adjustment
- 1.20 Extending arm
- 1.21 Video attachment

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e-mail: asistencia@auxilab.com

<http://www.auxilab.es>

Thank you for choosing a Colcoscope ZUZI. Its high mechanical quality will allow you to make use of this equipment for a long time.

Nevertheless, to get the best results and a longer life out of this equipment we recommend to read carefully this manual and follow the enclosed instructions.

In order to use this equipment, please take in mind the following points:

- ♦ This manual is inseparable from the Colcoscope ZUZI and it must be available to all the users of this equipment.
- ♦ Any doubts can be clarified by getting in touch with the ZUZI technical service.
- ♦ As it is a precision instrument, it should be handled carefully avoiding sudden movements and bumps.
- ♦ Never dismantle the different pieces of the Coloscope to repair them yourself. In addition to losing the warranty, it could lead to a faulty performance of this equipment.
- ♦ Do not use fluids that can liberate steam or form explosive or inflammable mixtures.

ATTENTION!

NO EQUIPMENT WILL BE REPAIRED BY OUR TECHNICAL SERVICE UNLESS IT IS PROPERLY CLEAN AND DESINFECTED.

INDEX OF LANGUAGES

Spanish	2-8
English	9-15

INDEX OF CONTENTS

1. USES OF THE INSTRUMENT	10
2. DESCRIPTION	10
3. TECHNICAL SPECIFICATIONS	11
4. SETTING UP.....	12
5. INSTRUCTIONS OF USE	12
6. USEFUL ADVICE	13
7. MAINTENANCE.....	13
8. TROUBLESHOOTING	14
9. ACCESSORIES	14
ANNEX I: WARRANTY.....	15